ANWENDERINFORMATIONEN

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
info@infumed.de

User Information

Non-sterile high-pressure tubes and high-pressure manifolds

This information must be read carefully before use in order to ensure safe further processing and use of the products and to avoid damage to patients and users as well as property damage due to improper further processing and use!

- The non-sterile high-pressure tubes are used exclusively for the subsequent packaging and sterilization and subsequent single use as non-invasive links between contrast medium injector and patient access. They are used, after packaging and sterilization, for the purpose of transfer or injection of contrast media and sodium chloride solution under high pressure (up to 1233 psi / 85 bar) and are intended for short-term use.
- The packaged and sterilized high-pressure manifolds are used for filling or flushing the catheter placed in the patient with contrast medium and sodium chloride solution and for invasive blood pressure measurement. As a non-invasive link between the patient connection and the contrast medium syringe, they are also intended for short-term use after packaging and sterilization.
- Any use beyond this is considered to be improper. Never use the medical device for other purposes! Infumed GmbH guarantees the suitability of the high-pressure tubes and manifolds for the intended use. Other manufacturer specific requirements must also be observed.

Packing and sterilization instructions:

- Before use, the non-sterile high-pressure tubes and manifolds must be packed in a sterile barrier system in accordance with EN ISO 11607-1.

- It is the responsibility of the processor to package the products under appropriate, monitored conditions and to monitor the microbiological and particulate contamination and endotoxins in accordance with the intended use.
- The packaging and the packaging process must be validated according to EN ISO 11607-2, regularly monitored and suitable for the subsequent sterilization process.

The following ETO sterilization parameters have been validated for the sterilization of the non-sterile high-pressure tubes and manifolds in the vacuum process:

<u>Ethylene oxide:</u>	100 Vol. %
<u>Preconditioning:</u>	
Duration:	10 - 48 hours
Temperature:	45°C - 53°C
rel. Humidity:	45% - 80%
<u>Sterilization:</u>	
Initial vacuum:	100 + / - 3 mbar
Initial vacuum reached:	8 - 60 minutes
Chamber leak rate:	max. 18 mbar in 6 minutes
Pressure increase when water vapor is input:	20 - 25 mbar
Duration of water vapor input:	2 - 30 minutes
Final pressure after	
ETO addition:	450 - 750 mbar
ETO dosage:	6.7 - 9.9 kg
Duration of ETO addition:	5 - 60 minutes
Chamber temperature:	45°C - 53°C
Exposure time:	4 - 5 hours
Rinsing cycles	3 x 200 + / - 3 bis 700 + / - 3 mbar 2 x 150 + / - 3 bis 850 + / - 3 mbar
Duration of rinsing cycles	45 - 165 minutes
<u>Postconditioning:</u>	
Degassing time:	22 - 54 hours
Temperature:	45°C - 53°C
Air change:	>7 times / h

- Under the above mentioned sterilization conditions, it has been validated and proven that the products maintain their sterility. It is the responsibility of the processor to ensure that the sterilization achieves the desired result under the selected process-specific conditions with the selected parameters and materials. This requires validation of the sterilization process and regular process monitoring.
- During packaging and sterilization, care must be taken to ensure that the products are not compressed or kinked.
- It is the responsibility of the processor to ensure, through the appropriate and validated packaging and sterilization process, that sterility is maintained under the storage and transport conditions intended by the processor until the sterile packaging is damaged or opened.
- Labelling: The processor is responsible for labelling the packaged and sterilized product. He is also responsible for providing instructions for use, warnings and precautions. The sterile packaging must be marked in accordance with the applicable regulatory requirements and the applicable standards. In particular, the specifications of EN ISO 15223-1, EN 1041, EN ISO 8536-4, -9, -10, -12 must be observed.
- The processor must supplement these instructions on use and application as well as warnings and cautions with his own instructions on packaging, sterilization and transport and pass them on to the user.
- The medical device may only be used by suitably qualified medical personnel, who are familiar with any complications that may arise. Note: Complications can occur at any time during use!
- The high-pressure tubes and manifolds may only be used in combination with commercially available consumables with Luer-Lock connections and with contrast medium injectors, whereby the specified maximum pressure of the respective product may not be exceeded.
- Remove the protective caps only immediately before use to ensure the sterility of the medical device until it is used.
- Indications, contraindications and use of the product must always be based on current medical standard. Guidelines and recommendations of the respective medical societies must be taken into account accordingly.
- To avoid contamination of the medical device, observe the general principles of hygiene according to the RKI (Robert Koch Institute) and the special hygiene requirements of your institution.
- The product is packaged in a protective manner and must be protected from moisture and direct sunlight and stored at a temperature between 5°C and 30°C until it is used.
- If stored properly, the product can be used up to the specified expiry date (see packaging).

- The sterile medical device is only intended for single use on one patient. It must not be processed, cleaned, disinfected or resterilized after use. In the event of reprocessing, the medical device would no longer be usable and would pose a potential risk to patients and users. Reprocessing can change the chemical, biological, physical or microbiological properties, which can result in a loss of functionality of the medical device, therefore the product must be discarded after single use. Due to the surface texture of the individual components of the medical product, successful cleaning after use on the patient cannot be guaranteed. A possible contamination is therefore unavoidable when used again. Resterilization can negatively affect the surface and material properties of the materials used. As a result, the original functionality and the product-specific integrity of the medical device can no longer be guaranteed.

- If the medical device itself is damaged, it must no longer be used and must be disposed of.

- We would like to point out that sterilization, storage and transport can cause individual screwed connections and protective caps to loosen. Before use, check all screwed connections / protective caps and tighten these connections if necessary.

- Make sure that combination products, that do not belong to the medical device, are compatible and follow their instructions for use.

- In the case of medical products with built-in check valves, it should be noted that the flow of the medium is only possible in one direction.

- Before use, it must be ensured that the medical device is free of air bubbles. An air embolism can lead to death or serious injury.

- The transparent plastic components used in our medical devices offer the advantage of visual control of the medium. When wetted with alcohol-based disinfectants or lipid-based solutions, stress cracks can occasionally occur on components made of plastic (e.g. made of polycarbonate or acrylonitrile-butadiene-styrene). Always avoid contact of disinfectant with the medical device.

- The medical device must be disposed of immediately after it has been used on a single patient. The relevant public and in-house guidelines must be observed for the disposal of medical products. After the medical device has been used, it can be potentially infectious and therefore a biological hazard may be present. Please observe your hygiene regulations so that contamination can be ruled out.

- Note: All serious incidents related to the device must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is established.

Please note the following symbols on the product labels and in the user information:

Symbol	Meaning
	Reference number
	Batch code
	Date of manufacture (year-month)
	Use-by date (year-month)
	Non-sterile
	Keep away from sunlight!
	Keep dry!
	Do not re-use!
	temperature limit of storage min. and max.

Symbol	Meaning
	Quantity per carton
	Consult instructions for use!
	Caution!
	Manufacturer
	CE marking with identification number (indicates that the medical device meets the requirements of the European Medical Device Directive 93/42/EWG)

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Infumed Med. Schlauchsysteme und Filtertechnik GmbH