DE ANWENDERINFORMATIONEN

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EN USER INFORMATION

High-pressure connectors and manifolds

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Ausgabe 3_11/2019

This information must be read carefully prior to use to ensure that the products are used safely and to prevent damage caused by misuse with patients and users, as well as damage to materials and equipment.

- These high-pressure connectors are non-invasive links between the contrast agent injector and the patient delivery point. They are used to percolate or inject contrast agents and solutions of sodium chloride under high pressure (up to 1233 psi / 85 bar) and are intended for short-term application.
- High-pressure manifolds are used to fill or to flush the catheter placed in the patient with a contrast agent and a solution of sodium chloride, and for invasive blood pressure measurement. They are also used for short-term non-invasive links between patient connection and contrast agent injector.
- Any other form of use is considered contrary to its intended use. Never use medical products for any other than their intended purpose! Infumed GmbH guarantees that these high-pressure connectors and manifolds are suitable for their intended use. Also pay attention to any other manufacturer-specific requirements.
- This medical product may only be used by appropriately qualified medical personnel familiar with the possible complications that can occur. Note: Complications can occur at any time during the use of this product.
- These high-pressure connectors and manifolds may only be used in combination with commercially available consumables with Luer-Lock connections and contrast agent injectors. Never exceed the stipulated maximum pressure for each product.
- Do not remove the protective caps until immediately before use to ensure the sterility of the medical product until it is used.
- Indications, contra-indications and use of the product must always be in accordance with the applicable medical standard. Always pay attention to the guidelines and recommendations of the relevant medical specialist companies.
- To prevent contamination of this medical product, pay attention to the general principles of hygiene as defined by RKI (Robert Koch Institute) and the particular hygiene requirements of your company.
- The product is packaged in a sterile condition (ethylene oxide sterilization) and must be protected from humidity and direct sunlight until it is used and must then be stored at a temperature of between 5°C and 30°C.
- When stored properly, the product can be used up to its stated use-by date (see packaging).
- The product is only sterile if the sterile packaging is undamaged and unopened and the specified use-by date has not expired. If the sterile packaging was damaged or accidentally opened before use, or if the use-by date has expired, the product must be disposed of.
- This sterile medical product is only intended for one-time use on one patient. It must not be processed, cleaned, disinfected or resterilized. If reprocessed, the medical product cannot be reused since it would constitute a potential hazard to patients and users. Reprocessing may alter the chemical, biological, physical or microbiological properties, causing the medical product to lose its function. Due

to the surface properties of individual components of the medical product, successful cleaning cannot be guaranteed after use on a patient. Consequently, if reused, contamination is an unavoidable possibility. Reesterilization can adversely affect the surface and material properties of the materials used. The original functionality and the product-specific integrity of the medical product can then no longer be guaranteed.

- If the medical product is damaged, it must not be reused and needs to be disposed of.
- We point out that sterilization, storage and transport can cause individual screwed connections and protective caps to become loose. Before use, check all screwed connections / protective caps and, if necessary, retighten these connections.
- Ensure that combination products not belonging to the medical product are compatible and pay attention to their instructions for use.
- For medical products fitted with non-return valves, you have to make sure that the medium can only flow in one direction.
- Before use, ensure that the medical product is free of air bubbles. An air embolism can cause death and can cause serious injuries.
- The transparent plastic components used in our medical products enable the medium to be checked visually. When wetted with disinfectants containing alcohol or solutions that contain lipids, isolated stress cracks can occur in plastic components (e.g. made of polycarbonate or acrylonitrile-butadiene-styrene). Always avoid contact between disinfectants and the medical product.
- This medical product must be disposed of after one-time use on one patient. Comply with applicable statutory and in-house directives when disposing of medical products. A bio-hazard can exist following the use of this medical product. Please pay attention to your hygiene regulations to eliminate any possibility of contamination.

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und Filtertechnik GmbH**